


DISPOSABLE

Endoscopy Instruments

Contributing to safety Disposable Instruments

Infection control in endoscopy has been a growing concern. Possible endoscopic transmission of infection has recently come under greater scrutiny. Single-use endoscopy instruments are attracting greater attention. Fujinon disposable Endoscopy Instruments offer high performance at low cost.



**Standard Fenestrated 2.8**


The hole in the center of each cup enables the sampling of more tissue while minimizing surrounding tissue damage.

Code	Needle	Sheath	Channel	Working Length	Qty
BF2316DF	-	-	φ2.8mm or larger	1600mm	20pieces/box
BF2318DF	-	-	φ2.8mm or larger	1800mm	

**Standard Fenestrated with Sheath 2.8**


The sheath reduces the scattering of mucus. The hole in the center of each cup enables the sampling of more tissue while minimizing surrounding tissue damage.

Code	Needle	Sheath	Channel	Working Length	Qty
BF2516DF1	-	yes	φ2.8mm or larger	1600mm	20pieces/box
BF2523DF1	-	yes	φ2.8mm or larger	2300mm	

**Standard Fenestrated Sheath with Needle 2.8**

The needle is set in the center of cup, preventing the cup from sliding.

Code	Needle	Sheath	Channel	Working Length	Qty
BF2516DN1	yes	yes	φ2.8mm or larger	1600mm	20pieces/box
BF2523DN1	yes	yes	φ2.8mm or larger	2300mm	

**Standard Fenestrated 2.0**

The hole incorporated in the center of each cup minimizes surrounding tissue damage.

Code	Needle	Sheath	Channel	Working Length	Qty
BF1716DF	-	-	φ2.0mm or larger	1600mm	20pieces/box


**Standard Fenestrated for Endoscopes**

These forceps are designed for small bowel biopsy. The hole incorporated in the center of each cup enables the sampling of more tissue while minimizing surrounding tissue damage. Sheath (optional) reduces the scattering of mucus.

Code	Needle	Sheath	Channel	Working Length	Qty
BF2525DF1	-	yes	φ2.8mm or larger	2500mm	5pieces/box
BF1725DF	-	-	φ2.0mm or larger	2500mm	

Washing Brush

Brushes for cleaning inside of the channel.

Washing Brush**Channel Washing Brush**

To facilitate channel washing for wide-scope forceps.

Code	Channel	Working Length	Qty
WB3225DS	φ2.0mm ~ φ2.8mm	2500mm	20pieces/box
WB4325DS	φ3.2mm ~ φ4.2mm	2500mm	

FUJIFILM GROUP

FUJINON CORPORATION

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN

TEL : 81-48-668-2153 FAX : 81-48-668-1570

E-mail : sales@msv.fujinon.co.jp

URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.

TEL : 1-800-490-0661 FAX : 1-973-633-8818

E-mail : med.info@fujinon.com

URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILlich, GERMANY

TEL : 49-2154-924-0 FAX : 49-2154-924-290

E-mail : Fujinon@fujinon.de

URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK, SINGAPORE, 159552

TEL : 65-6276-4988 FAX : 65-6276-6911

E-mail : fujinon@fujinon.com.sg

URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA

TEL : 61-2-9748-2744 FAX : 61-2-9748-2428

E-mail : sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO., LTD.

FLAT B, 6/FL, SHI MEI BLD, 79 RI JING ROAD,

WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIPCODE: 200131 CHINA

TEL : 86-21-5424-9262(Office) TEL : 86-21-5866-2151(Service)

FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com