

FUJINON Video Bronchoscope

New

EB-530H/S/T

Fujinon's most advanced 500 series Video Bronchoscopes are designed for physician productivity and patient care.



FUJINON
FUJIFILM



FUJINON Video Bronchoscope

Excellent image by Super CCD technology
 Ideal Flexibility of New Insertion Shaft
 Hygienic Disposable Button
 Compatible with FICE (Flexible spectral-Imaging Color Enhancement)

EB-530S

Standard Type

Well balanced bronchoscope suitable for standard examination.

EB-530T

Treatment Type







2.8mm channel for various accessories for treatment

EB-530H

Super image Type

High resolution with big brilliant image especially matching to FICE observation.

Specifications

Model	EB-530S	EB-530T	EB-530H	EB-530S	EB-530T	EB-530H
Field of view	120°	120°	140°			
Observation range	3-100mm	3-100mm	3-100mm			
Distal end diameter	4.9mm	5.9mm	5.4mm			
Flexible portion diameter	4.9mm	5.9mm	5.1mm			
Channel diameter	2.0mm	2.8mm	2.0mm			
Bending capability: UP/DOWN	180°/130°	180°/130°	180°/130°			
Working length	600mm	600mm	600mm			
Total length	870mm	870mm	870mm			

FUJIFILM Group

FUJINON Corporation

1-324 Uetake, Kita-ku, Saitama City, Saitama 331-9624, Japan
 TEL: 81-48-668-2153 FAX: 81-48-668-1570
 URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 High Point Drive Wayne, NJ. 07470, U.S.A.
 TEL: 1-800-490-0661 FAX: 1-973-633-8818
 E-mail: med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

Halskestrasse 4, 47877 Willich, Germany
 TEL: 49-2154-924-0 FAX: 49-2154-924-290
 E-mail: fujinon@fujinon.de URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE.LTD.

Blk 211 Henderson Road, #10-04 Henderson Industrial Park, Singapore 159552
 TEL: 65-6276-4988 FAX: 65-6276-6911
 E-mail: fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.

Unit 18, 52 Holker Street, Silverwater N.S.W 2128, Australia
 TEL: 61-2-9748-2744 FAX: 61-2-9748-2428
 E-mail: sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO., LTD.

Flat K, 23/F, 985 Dong Fang Road, Pudong District, Shanghai 200122, China
 TEL: 86-21-5424-9262 FAX: 86-21-5424-9264
 E-mail: fing2001@hotmail.com