


FUJINON Video Colonoscope

New

EC-530
FM/FI/FL



FUJINON
FUJIFILM



FUJINON Video Colonoscope


EC-530FM EC-530FI EC-530FL

For the Lower G.I. Tract – Standard Type

These super wide-angle standard colonoscopes offer a large 3.8mm working channel inside a slim 12.8mm outside diameter. An ultra-wide 140 degree field of view enhances image quality from the center to the edge. These colonoscopes also offer a wider observation range from 3-100mm. In addition, an integrated forward water jet allows for lavage in clinical situations.

Specifications

Model	EC-530FM/FI/FL
Field of view	140°
Observation range	3-100mm
Distal end diameter	12.8mm
Flexible portion diameter	12.8mm
Bending capability: UP/DOWN	180°/180°
LEFT/RIGHT	160°/160°
Working length	1,330/1,520/1,690mm
Total length	1,630/1,820/1,990mm
Forceps channel diameter	3.8mm



FUJIFILM Group
FUJINON Corporation
 1-324 Uetake, Kita-ku, Saitama City, Saitama 331-9624, Japan
 TEL: 81-48-668-2153 FAX: 81-48-668-1570
 URL <http://www.fujinon.co.jp/>

FUJINON INC.
 10 High Point Drive Wayne, NJ. 07470, U.S.A.
 TEL: 1-800-490-0661 FAX: 1-973-633-8818
 E-mail: med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH
 Halskestrasse 4, 47877 Willich, Germany
 TEL: 49-2154-924-0 FAX: 49-2154-924-290
 E-mail: fujinon@fujinon.de URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE.LTD.
 Blk 211 Henderson Road, #10-04 Henderson Industrial Park, Singapore 159552
 TEL: 65-6276-4988 FAX: 65-6276-6911
 E-mail: fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.
 Unit 18, 52 Holker Street, Silverwater N.S.W 2128, Australia
 TEL: 61-2-9748-2744 FAX: 61-2-9748-2428
 E-mail: sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO., LTD.
 Flat K, 23/F, 985 Dong Fang Road, Pudong District, Shanghai 200122, China
 TEL: 86-21-5424-9262 FAX: 86-21-5424-9264
 E-mail: fing2001@hotmail.com