

# FUJINON Video Gastroscope

**New**

---

## **EG-530FP**

---



**FUJINON**  
FUJIFILM



# FUJINON Video Gastroscope

## EG-530FP

### For the Upper G.I. Tract - Slim Type

Fujinon's EG-530FP incorporates the most recent advancements in the miniaturization of CCD chips which offers a slim 8.5mm outside diameter and a 2.8mm working channel. A slim 8.5mm outside diameter and a 2.8mm working channel allow for easy insertion and standard accessory passage.

## Specifications

Field of view	140°	
Observation range	3-100mm	
Distal end diameter	8.5mm	
Flexible portion diameter	8.5mm	
Bending capability: UP/DOWN	210°/90°	
LEFT/RIGHT	100°/100°	
Working length	1,100mm	
Total length	1,400mm	
Forceps channel diameter	2.8mm	

FUJIFILM Group

### FUJINON Corporation

1-324 Uetake, Kita-ku, Saitama City, Saitama 331-9624, Japan  
TEL: 81-48-668-2153 FAX: 81-48-668-1570  
E-mail: [med.info@fujinon.com](mailto:med.info@fujinon.com) URL <http://www.fujinon.co.jp/>

### FUJINON INC.

10 High Point Drive Wayne, NJ. 07470, U.S.A.  
TEL: 1-800-490-0661 FAX: 1-973-633-8818  
E-mail: [med.info@fujinon.com](mailto:med.info@fujinon.com) URL <http://www.fujinonendoscopy.com/>

### FUJINON (EUROPE) GmbH

Halskestrasse 4, 47877 Willich, Germany  
TEL: 49-2154-924-0 FAX: 49-2154-924-290  
E-mail: [fujinon@fujinon.de](mailto:fujinon@fujinon.de) URL <http://www.fujinon.de/>

### FUJINON SINGAPORE PTE.LTD.

Blk 211 Henderson Road, #10-04 Henderson Industrial Park, Singapore 159552  
TEL: 65-6276-4988 FAX: 65-6276-6911  
E-mail: [fujinon@fujinon.com.sg](mailto:fujinon@fujinon.com.sg) URL <http://www.fujinon.com.sg/>

### FUJINON AUSTRALIA PTY. LTD.

Unit 18, 52 Holker Street, Silverwater N.S.W 2128, Australia  
TEL: 61-2-9748-2744 FAX: 61-2-9748-2428  
E-mail: [sales@fujinon.com.au](mailto:sales@fujinon.com.au)

### SHANGHAI FUJINON MEDICAL EQUIPMENT CO., LTD.

Flat K, 23/F, 985 Dong Fang Road, Pudong District, Shanghai 200122, China  
TEL: 86-21-5424-9262 FAX: 86-21-5424-9264  
E-mail: [fing2001@hotmail.com](mailto:fing2001@hotmail.com)