

FUJINON Video Gastroscope

New

EG-530N
EG-530NP



FUJINON
FUJIFILM



FUJINON Video Gastroscope



EG-530N EG-530NP

For the Upper G.I. Tract-Super Slim scope

Fujinon's new super slim trans-nasal endoscope provides for a higher level of trans-nasal endoscopy.

An impressive 4.9mm (EG-530NP) and 5.9mm(EG-530N) distal tip diameter, combined with a 2.0mm channel, and 2 light guides for better illumination allow for improved insertion capability.

Specifications

	EG-530N	EG-530NP	EG-530N	EG-530NP
Viewing direction	Forward	Forward		
Observation range	3-100mm	3-100mm		
Field of view	120°	120°		
Distal end diameter	5.9mm	4.9mm		
Flexible portion diameter	5.9mm	5.1mm		
Bending capability: UP	210°	210°		
DOWN	90°	120°		
LEFT/RIGHT	100°			
Forceps channel diameter	2.0mm	2.0mm		
Working length	1,100mm	1,100mm		
Total length	1,400mm	1,400mm		

FUJIFILM Group

FUJINON Corporation

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN
TEL : 81-48-668-2153 FAX : 81-48-668-1570
E-mail:sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.
TEL : 1-800-490-0661 FAX : 1-973-633-8818
E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILlich, GERMANY
TEL : 49-2154-924-0 FAX : 49-2154-924-290
E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,
SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911
E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA
TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD,79 RI JING ROAD,
WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA
TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)
FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com