

FUJINON Digital Video Processor

New

EPX-2500



FUJINON
FUJIFILM

FUJINON Digital Video Processor

EPX-2500


Digital Processor and 150W Xenon Light Source Combination Unit

Fujinon introduces a higher standard for electronic video endoscopy with the newly designed EPX-2500. Setting a new criteria for diagnostic and therapeutic endoscopy, the EPX-2500 provides physicians with Fujinon's world-class imaging and the exceptional ergonomics of the newly designed 500 series endoscopes.

Specifications

Digital Output	DVI (Digital Visual Interface) : 1024x768p
Analog Output	RGB(2) : SDTV (NTSC / PAL) Y/C(2) : SDTV (NTSC / PAL) Composite : SDTV (NTSC / PAL)
Color Adjustment	Black, Red, Green, Blue, R-Hue, Chroma; 9 settings
Detail	Hi, LO; 9 settings
Contrast (Gamma)	9 settings
BLV	Hi, Mid, Lo, Off
Picture in Picture	On, Off; Size: 1/4, 1/3
Auto Gain Control	Off, +3db, +6db
Iris	Average/Peak
Zoom	Electric zoom : x1.0 - x2.0; 0.05steps
Lamp rated value	Main lamp : 11.7V 150W Xenon lamp Emergency lamp : 12V 75W Halogen lamp
Brightness control	9 settings
Lamp cooling method	Forced air cooling
Air supply pump	Hi, Low, Off
Power	120V 60Hz 2.7A / 230V 50Hz 1.4A
Dimensions (W x H x D)	375 x 190 x 495mm
Weight	17.0kg

Function



FUJIFILM Group

FUJINON Corporation

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN
TEL : 81-48-668-2153 FAX : 81-48-668-1570
E-mail:sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.
TEL : 1-800-490-0661 FAX : 1-973-633-8818
E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY
TEL : 49-2154-924-0 FAX : 49-2154-924-290
E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,
SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911
E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA
TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD,79 RI JING ROAD,
WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA
TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)
FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com