

FUJINON

Video Cysto Endoscopes



FUJINON
 **FUJIFILM**

FUJINON VIDEO CYSTO ENDOSCOPES


EY-470S/B & EY-470S

Fujinon is proud to announce its NEW Video Cystoscopes, the EY-470S/B and the EY-470S -- a giant leap forward in surgical technology. Through a marriage of proprietary video chip technology and an industry-leading flexible scope design, Fujinon has created flexible video cystoscopes with a higher resolution image than ever before possible. Images are large, clear and consistently bright. These NEW flexible video cystoscopes are ideal for the office, clinic or hospital environment whether performing diagnosis or, in the case of the EY-470S, treatment with a new suction feature that can be controlled by one hand to maintain a clear view.



- High resolution image resulting in better observation through a high-density CCD chip.
- Large instrument channel for better flow and wider range of endoscopic instruments.
- Suction capability (available with the EY-470S) that can be controlled to maintain clear view.
- HF compatibility enables Urologists to perform therapeutic endoscopy.
- A tapered distal end greatly improves insertion capability.

	EY-470S/B	EY-470S
Viewing direction	Forward	Forward
Observation range	3 ~ 100mm	3 ~ 100mm
Field of view	120°	120°
Distal end diameter	5.1mm	5.1mm
Flexible portion diameter	5.1mm	5.1mm
Bending capability: UP / DOWN	210° / 130°	210° / 130°
Forceps channel diameter	2.0mm	2.0mm
Working length	380mm	380mm
Total length	650mm	650mm
Suction	Not available	Available



FUJIFILM GROUP

FUJINON CORPORATION

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN
 TEL : 81-48-668-2153 FAX : 81-48-668-1570
 E-mail:sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.
 TEL : 1-800-490-0661 FAX : 1-973-633-8818
 E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY
 TEL : 49-2154-924-0 FAX : 49-2154-924-290
 E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,
 SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911
 E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA
 TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD,79 RI JING ROAD,
 WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA
 TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)
 FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com