

# FUJINON Fiber Choledocho Endoscopes



**FUJINON**  
 **FUJIFILM**

## FUJINON FIBER CHOLEDOCHO ENDOSCOPES

# FO-120S

Fujinon is renowned for its versatility and advanced technology in all areas of endoscopy. In order to meet the diverse needs of today's healthcare providers, Fujinon has developed a new specialty scope in the field of General Surgery, the FO-120S Fiber Choledochoscope.



- The FO-120S is a truly unique endoscope that incorporates a large 2.2mm working channel within a slim 5.3mm insertion section. The slim size increases maneuverability and allows for easy introduction even through areas with severe strictures.
- Thanks to an ultra-small, high density fiber bundle, the FO-120S enables you to see the highest resolution images. This advancement results in outstandingly crisp images.
- In order for the endoscope control section to be held more comfortably, a newly designed lightweight control section has been designed. This innovation also lessens hand fatigue during the operation.

	FO-120S	
Viewing direction	Forward	
Observation range	1 ~ 50mm	
Field of view	120°	
Distal end diameter	4.8mm	
Flexible portion diameter	5.3mm	
Bending capability : UP	180°	
DOWN	130°	
Forceps channel diameter	2.2mm	
Working length	450mm	
Total length	750mm	

FUJIFILM GROUP

### FUJINON CORPORATION

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN  
 TEL : 81-48-668-2153 FAX : 81-48-668-1570  
 E-mail:sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

### FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.  
 TEL : 1-800-490-0661 FAX : 1-973-633-8818  
 E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

### FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY  
 TEL : 49-2154-924-0 FAX : 49-2154-924-290  
 E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

### FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,  
 SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911  
 E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

### FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA  
 TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

### SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD,79 RI JING ROAD,  
 WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA  
 TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)  
 FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com