

# FUJINON Light Sources



## HALOGEN LIGHT SOURCE PS2-HS

A 150 Watt halogen light source, excellent for diagnostic use. A spare lamp is built-in so that should the primary lamp fail, during a procedure, no interruption of the examination will occur. An air pump, for insufflation, is incorporated, and an on/off switch is provided for procedures which do not require it. Light intensity is manually controlled.



## HALOGEN LIGHT SOURCE PS2-HP

This sturdy, reliable light source is designed for its portability. A single 150 Watt halogen lamp, whose light intensity is manually controlled, is used. An air pump, for insufflation, is incorporated, and an on/off switch is provided for procedures which do not require it.



## COMPACT HALOGEN LIGHT SOURCE PS2-HN

Lightweight and compact, this 150 Watt halogen light source is designed for diagnostic procedures such as Bronchoscopy, Hysteroscopy and Nasopharyngoscopy, where no air pump is required. It has a built-in camera drive mechanism, and a manual light intensity control.

	PS2-HS	PS2-HP	PS2-HN
Power supply	AC117V / 230V		
Frequency	50 / 60Hz		
Bulb used	15V, 150W halogen bulb, gold mirror		
Spare bulb	Incorporated	-	
Air pump	○		-
Dimensions (W) x (D) x (H)	340 x 420 x 145mm	195 x 350 x 215 mm	195 x 350 x 215 mm
Weight	12 kg	8.6 kg	7.9 kg

# FUJINON


**FUJIFILM**

FUJIFILM GROUP

### FUJINON CORPORATION

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN  
TEL : 81-48-668-2153 FAX : 81-48-668-1570  
E-mail:sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

### FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.  
TEL : 1-800-490-0661 FAX : 1-973-633-8818  
E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

### FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY  
TEL : 49-2154-924-0 FAX : 49-2154-924-290  
E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

### FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,  
SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911  
E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

### FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA  
TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

### SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD,79 RI JING ROAD,  
WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA  
TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)  
FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com